



Black Cube SBS Rack for Pipetting Tasks

Description

Black Cube is a stable SBS rack and suitable for all suitable for all standard formats. The Black Cube is filled with water and placed in a freezer before use. placed in a freezer before use.

The material mix of sturdy plastic and a metal base provides an ideal receptacle for PCR plates in automated liquid handling systems.

A modern and versatile rack suitable for many applications.

To prepare, freeze the PCR Cooler at -20°C for at least 2 hours. Turn the bottom side upwards when doing so.

Filling quantity: 350 ml
Optional: PCR-Cooler plate made of aluminium (blue).

Base plate removable for lower stands.



SBS-Rack Black Cube



SBS-Rack Black Cube with PCR-Cooler (optional)



Black-Cube with removable base

Technical Data

Format	SBS standard
Water tank filling quantity	350 ml
Dimensions exterior (L x W x H)	135 x 94 x 110 mm
Dimensions of Mikroplate area (L x W x H)	127 x 85 x 15 mm
Dimensions removable base (L x W x H)	135 x 94 x 20 mm